

High resolution CESI-MS analysis of APTS-labeled N-glycans of biopharmaceutical interest

Bryan Fonslow¹, Boglarka Donczo², Marton Szigeti², Andras Guttman^{1,2}
¹SCIEX, Brea, CA and ²Horvath Csaba Laboratory of Bioseparation Sciences, University of Debrecen, Hungary

ABSTRACT

Glycosylation plays an important role in cellular interactions, protein folding, and monoclonal antibody (mAb) stability. For characterization, N-linked glycosylation on mAbs is routinely analyzed by CE-LIF after endoglycosidase based glycan release and fluorophore labeling. CE separation provides structural resolution and migration time-based identification of glycans, but coupling with MS detection offers additional structural information. The integration of CE and ESI into a single dynamic process (CESI) provides the capability of performing CE separation and MS ionization with ultra-low flow rates, resulting in reduced ion suppression and improved sensitivity. CESI-MS has been optimized and evaluated for APTS-labeled mAb and formalin-fixed paraffin-embedded glycan analysis.

INTRODUCTION

Since introduced as a biological tissue sample fixative by F. Blum in 1893, formalin fixation has been still an almost exclusively used method for tissue preservation. Fixation using this reagent makes clinically relevant samples stable that can be stored at ambient condition for decades in depositories. Combined with paraffin embedment, formalin fixed tissues become FFPE specimens (formalin fixed paraffin embedded) that are amenable for long-term storage and can be universally used in histopathological laboratories. The clinicopathological significance of formalin-fixed paraffin-embedded (FFPE) archival specimens has led many researchers to utilize and annotate these “unmasked treasures” beyond their main histopathological utilization even at the molecular level. Several clinico-molecular assays for oncological prognosis and treatment decision have been developed based on FFPE tissue analysis. Recently capillary electrophoresis has been introduced as a method to probe the N-glycosylation pattern of these precious samples.¹

MATERIALS AND METHODS

Sample Preparation: A maltooligosaccharide ladder, IgG1 glycan standards, and released formalin-fixed paraffin-embedded (FFPE) glycans were labeled with APTS and cleaned-up. For deparaffinization, the FFPE samples were washed twice by xylene for 20 min and centrifuged at 5000 × g for 10 min. This was followed by ethanol wash, two times for 20 min and centrifugation at 5000 × g for 10 min. For protein extraction, all samples were first homogenized in a 1 mL Tapered Tissue Grinder (VWR, Budapest, Hungary) and solubilized with 100 µL of RIPA buffer (25 mM Tris-HCl, pH 7.6, 150 mM NaCl, 1% NP-40, 1% sodium deoxycholate, 0.1% SDS) [8] containing 10 µL of 50 mM DTT. The solubilization mixture was first incubated at 100° C for 20 min, then 65° C for 60 min. After the solubilization step, 10 µL of 50 mM iodoacetamide was added and incubated at 37° C for 30 min in dark. For buffer exchange, 10 kDa spin filters were used at 5000 × g for 10 min. To release the N-glycans, 49 µL of 20 mM NaHCO₃ buffer (pH 7.0) and 1 µL PNGase F Ultra enzyme (2.5 mU, Prozyme, Hayward, CA) were added to the sample on the filters and incubated for 1 h at 50° C. The released N-glycans were centrifuged through the spin filters at 5000 × g for 10 min and dried in a centrifugal vacuum evaporator. 6 µL of 20 mM APTS and 2 µL of 1 M NaCNBH₃ (in THF) was added to the dry pellet and incubated at 37° C overnight. The labeled samples were magnetic bead purified (CleanSeq, Beckman Coulter, Indianapolis, IN) and immediately used or stored at -20° C for later analysis.

CESI 8000 Plus Stand-Alone CE Mode Conditions: The separations were monitored by laser induced fluorescence (LIF) detection using a 488 nm solid state laser with a 520 nm emission filter. 50 cm effective (60 cm total) length with 50 µm id NCHO capillaries were used filled with NCHO gel buffer system (both from SCIEX) in reversed polarity mode (cathode at the injection side). The samples were injected by 1 psi pressure for 5 s.

CESI 8000 Plus MS Mode Conditions: CESI experiments were carried out with a SCIEX CESI 8000 Plus system equipped with a temperature controlled autosampler and a power supply with the ability to deliver up to 30 kV. CESI-MS analyses were performed using a SCIEX CESI 8000 Plus system. CESI-MS experiments also employed a TripleTOF® 6600 system. CESI-MS separations were achieved using an OptiMS bare fused-silica capillary cartridge by application of -30 kV and low pressures (2 for 7.5 min and 5 psi for 10 min). Ammonium acetate (7.5 mM, pH 4.5) with 10% isopropanol (IPA) was used as the background electrolyte (BGE). Samples were dissolved in 10% BGE to 10 µL and loaded into nanoVials (P/N 5043467) for ~15 nL pressure injections.

MS/MS Conditions: A SCIEX TripleTOF® 6600 system with a NanoSpray® III source and CESI adapter were used. Negative ESI was performed at -1750 V scanning from 250 - 2000 m/z with 250 ms (MS).

Data Analysis: High resolution MS and MS/MS spectra were analyzed using SCIEX PeakView®. Glucose unit (GU) values were calculated by the GUcal software (<http://www.gucal.hu>) and preliminary assignment of the glycan structures of interest were aided by the Glycobase ver 3.0 from NIBRT (Dublin, Ireland).




Figure 1. CESI 8000 Plus High Performance Separation-ESI Module coupled to a TripleTOF® 6600 System. The core of the CESI-MS technology is an etched porous capillary tip that allows for electrospray ionization without dilution of the BGE and sample by a make-up liquid. Instead a conductive liquid is used to apply the ESI voltage through the porous tip.

RESULTS

The transition from current CE-LIF APTS-glycan analysis methods to CESI-MS requires the use of MS-compatible buffers. Separations of an APTS-labeled carbohydrate ladder and mAb glycan standards were analyzed to evaluate the general capabilities of CESI-MS for the analysis of APTS-labeled released glycans from proteins of interest. Separation and detection of the carbohydrate ladder and mAb glycan standards illustrated efficient separations (~6 sec peak widths) and sensitive detection in negative ESI mode.




Figure 2. CE-LIF analysis of (A) an APTS-labeled carbohydrate ladder and (B) a fucosylated series of mAb glycan standards.

Separations of mAb glycan standards using a MS-compatible ammonium acetate buffer (pH 4.5) were also compared to the industry-accepted CE-LIF analyses using the stand-alone CE functionality of the CESI 8000 Plus instrument. Very similar separation profiles were observed, indicating ease of method transfer and correlation.




Figure 3. CE-LIF analysis of mAb glycan standards (A) with non-MS compatible BGE and (B) with MS-compatible 7.5 mM ammonium acetate, pH 4.5 with 10% isopropanol.

The CESI-MS method was evaluated using real samples with additional sample preparation. Prior to CESI-MS experiments, the released glycan samples were first characterized with industry-accepted CE-LIF analyses using the stand-alone CE functionality of the CESI 8000 Plus instrument. FFPE glycans from different mouse tissues were analyzed and compared. Glycan sequencing was performed with a combination of glycolytic enzymes to determine the carbohydrate content of each glycan.




Figure 4. (A) CE-LIF profiles of the released and APTS labeled N-glycans from different mouse tissue samples. A: lung B: brain C: heart D: spleen E: liver F: kidney G: tumor tissue H: intestine. (B) Exoglycosidase array based sequencing of mouse lung N-glycans. Traces: sialidase (★), sialidase + fucosidase (◆★), sialidase + fucosidase + galactosidase (◆◆★) and sialidase + fucosidase + galactosidase + N-acetylglucosamine (◆◆◆★).

Analysis of FFPE mouse lung and brain glycans were then performed by CESI-MS. Similar numbers of glycan peaks were extracted from the MS data for each sample. Accurate molecular masses were calculated from m/z's and charge states to find potential carbohydrate compositions for each APTS-glycan electrophoretic peak.




Table 1. List of APTS-glycan masses measured between FFPE mouse lung and brain tissue samples. Potential glycan compositions were assigned based on accurate molecular weights for the majority of measured masses.

CONCLUSIONS

The combined CESI-MS and stand-alone CE functionality of the CESI 8000 Plus system provides powerful methods to characterize and quantify released glycans from biologics and FFPE samples. Both analyses deliver high resolution separations, structural characterization, and relative quantitation. When coupled with the TripleTOF® 6600 system, the high sensitivity CESI-MS analyses are achieved from small sample amounts (~10 ng) and volumes (1 – 5 µL) using the nanoVial. Further analyses by MS/MS would facilitate confirmation of glycan structures and correlation to CE-LIF sequencing results.

REFERENCES

- Donczo, B., Szigeti, M. Ostoros, G., Gacs, A., Tovari, J., Guttman, A. N-Glycosylation analysis of formalin fixed paraffin embedded samples by capillary electrophoresis, *Electrophoresis*. **2015**.

TRADEMARKS/LICENSEING

For Research Use Only. AB Sciex is doing business as SCIEX. Not for use in diagnostic procedures. The trademarks mentioned herein are the property of AB Sciex Pte. Ltd. or their respective owners. © 2016 AB Sciex.