

Complete characterization of biotherapeutic proteins by automated data processing on high resolution accurate mass spectrometry with SWATH® acquisition



Milla Neffling¹; Yang Shi²; Dominic Gostick²; Doug Simmons³

¹SCIEX, Warrington, UK; ²SCIEX, Framingham, MA, USA; ³SCIEX, Concord, ON, Canada

INTRODUCTION

Biotherapeutic proteins are a constantly growing market in the pharmaceutical industry. The rise of the biosimilars market has increased the number of laboratories requiring robust analytical capabilities. Confidence in data quality and reproducibility is vital for the successful characterisation of biotherapeutic proteins. Confidence is not only dependent on the data quality, but also the ease-of-use in which it is obtained to avoid mistakes and the need for reacquiring the data.

Biotherapeutic proteins carry more complexity than the traditional small molecule drugs. The heterogeneity assessments are vital to ensure product quality – drug efficacy and patient safety. Mass spectrometry techniques are a powerful tool to answer the critical questions: full amino acid sequence, post-translational modification (PTM) ratios, as well as possible process and storage derived modifications.

Mass spectrometry techniques are used to investigate the biotherapeutic proteins on the level of intact, subunit and peptide mapping. Single platform able to handle all levels of characterisation is desired. This approach has been taken in this study to characterise two commercially available monoclonal antibodies (mAbs; see figure 1.) Rituximab and Trastuzumab. The data acquisition has been performed with a quadrupole-time-of-flight (Q-TOF) instrument. Here we utilise for peptide mapping for the first time SWATH® acquisition, a data independent acquisition (DIA) technique, unique to SCIEX instrumentation. The SWATH technique overcomes one of the issues associated with the default data acquisition methods: peptides (precursor ions) being irreproducible picked for MS/MS analysis. With SWATH acquisition all peptides are selected for MS/MS analysis – without instrument (data-dependent) or prior user (targeted) precursor ion selection.

SWATH acquisition has been previously used in the proteomics world to quantitatively analyse samples with high complexity (1). For this purpose, the data processing has relied on an ion library. With biotherapeutic proteins, the expectation is to have limited number of high abundant proteins in the samples – in which the low abundant proteoforms need to be assessed with confidence. With this in mind, the data processing for BioPharmaView™ software 2.0 has taken a different approach – using no ion library, but matching the theoretical MS/MS fragments.




Fig. 1

Figure 1. Large numbers of biotherapeutic proteins belong to the class of monoclonal antibodies. These proteins carry two heavy and two light chains covalently bound by disulphide bonds. The usual characterization workflows include full intact analysis for the approximately 150kDa protein complex, subunit analysis achieved by cleaving the protein by the hinge region, and/or by reducing the disulphide bonds, and by peptide mapping after a digestion with a proteolytic enzyme.

MATERIALS AND METHODS

1) intact protein analysis; 2) subunit analysis; and 3) peptide mapping analysis

Sample Preparation:

mAb samples were prepared following standard procedures for sample preparation. Shortly:

1) mAb material was diluted with aqueous 0.1% Formic Acid solution to a concentration of 0.02µg/µl.

2) For the subunit analysis the material was digested with Promega IdeZ protease according to manufacturer's recommendation (37°C – 1h). For Light chain – Heavy chain analysis the sample was reduced with TCEP (SCIEX protein digestion kit)

3) The tryptic digestion used octyl-β-D-glucopyranoside (OGS) as a denaturant, TCEP as a reducing agent, and MMTs as the alkylating reagent. All reagents were from an in-house SCIEX trypsin digest reagent kit and were used according to the instructions in the kit.

LC-UV-MS Conditions:

Two mAbs were analysed using as instrumentation a Shimadzu analytical flow instrument coupled to a quadrupole-time-of-flight accurate mass instrument. The instrument was operated in TOF MS, IDA and SWATH modes.

The chromatographic conditions used were standard for the respective workflows, utilising 1) and 2) Agilent C-8 Poroshell column with 15 to 20 min gradient run time, and 3) Phenomenex peptide Aeris C-18 column with 60 min runtime. The UV was collected at 214nm and 280nm.

Data processing

The data analysis (1)+ 2) data reconstruction and 3) peptide mapping) was performed with BioPharmaView™ software 2.0. System performance on mass accuracy, S/N levels and quantitative information on post-translational modifications (including glycosylation) were used to assess the reproducibility of the analysis.

Figure 2. Data processing was achieved with BiopharmaView software 2.0 which utilizes projects: in the project you specify the assay information (sequence, disulphide bonds and modifications) and characterize a golden standard. This processing is then utilized to number of samples in a batch to achieve fast and easy comparison studies of biotherapeutic samples (batches, stability study time points, biosimilarity assessments). The results can then be reviewed and a report generated.




Fig. 2

RESULTS

The current platform performance was demonstrated on several levels of analysis:

1) Intact protein, 2) subunit analysis and by 3) peptide mapping. The workflows allow for different levels of information:

