

Carbonic anhydrase inhibitory potential of *Lagenaria siceraria* and identification of its bioactive compounds-An LC-MS/MS approach



Joydeb Chanda¹, Pulok K Mukherjee^{1*}, Rajarshi Biswas¹, Dipankar Malakar², Manoj Pillai²

¹School of Natural Product Studies, Department of Pharmaceutical Technology, Jadavpur University, Kolkata -700 032, India.; ²Sciex, 121 UdyogVihar, Gurgaon, Haryana, India

ABSTRACT

Lagenaria siceraria Stand is an important member of Cucurbitaceae family, widely consumed as a vegetable in daily food habits among Indians. The current study was aimed to identify the active chemical compounds of *L. siceraria* involved in carbonic anhydrase inhibition through bioassay and liquid chromatography-mass spectrometry guided analysis. The extraction of fruits of *L. siceraria* was carried out in methanol. Fractionation of the methanolic extract *L. siceraria* (LSME) was performed in three different solvents (hexane, ethyl acetate and water) successively. The carbonic anhydrase inhibitory activity was evaluated through its esterase inhibition assay. Liquid chromatography and mass spectroscopy (LC-MS/MS) analysis of the bioactive fraction was performed to identify the major bio-actives present in it. The carbonic anhydrase inhibition assay revealed that the LSAF (aqueous fraction) possess highest esterase inhibition activity. The IC₅₀ value of LSAF was found equipotent activity with acetazolamide. LC-MS/MS analysis of LSAF revealed the presence of phenolic compounds, confirmed by their MS/MS spectrum. This finding suggests that phenolic compounds of *L. siceraria* may be responsible for carbonic anhydrase inhibition activity.

INTRODUCTION

Lagenaria siceraria Stand popularly known as Bottle gourd is an important food plant of cucurbitacea family. It is very well conversed in Ayurveda and folk medicine for its some potential therapeutic purposes. The phyto-constituents found in the fruit includes cucurbitacin B, phenolic glycosides, flavonoids, flavon-C-glycoside such as isovitexin, isoorientin, saponarin, phenolic acids like caffeic acid, cinnamic acid, flavonoids, sterols like fucosterol, campesterol, cytotoxic polysaccharides (Gangwal et al., 2010, Jaiswal et al., 2014). The plant possesses several significant therapeutic potentials include diuretic, antioxidant, antihyperglycemic, antihyperlipidemic, cardioprotective activity, immunomodulatory, anthelmintic activity (Ghule et al., 2009). Several phenolic compounds (phenolic acids, flavonoids) are considered beneficial for their activity. Carbonic anhydrase is a metalloenzyme, is involved in several pathophysiological processes in human. The present work was undertaken to establish the carbonic anhydrase inhibitory property of the bioactive fraction of LS and identifying its major phytoconstituents by LC-MS/MS study. This strategy was found to be very useful for rapid identification of bioactive components in plant extract prior to isolation of the pure compounds.

MATERIALS AND METHODS

Carbonic anhydrase II from bovine erythrocytes (3848 W/A units/mg solid) (EC-232-576-6) (bCA II) and p-nitrophenyl acetate (p-NPA) was purchased from Sigma Aldrich, St. Louis, MO, USA. Acetazolamide IP (Batch No. AZM-V-P/131107) was procured as a gift sample from Mangalam Drugs and Organics Ltd., Mumbai, India. Other chemicals and solvents were purchased from Merck Pharmaceuticals, Mumbai.

In-vitro carbonic anhydrase activity:

In-vitro carbonic anhydrase inhibition assay was performed through esterase inhibition model (Verpoorte et al., 1967). The absorbance of the samples in each well was determined at 400 nm using a UV-visible spectroscopy (SpectraMax Plus, United States). The change of absorbance was observed due to the liberation of p-nitrophenol as the hydrolysis product of p-NPA. The assay procedure was carried out in triplicate. Acetazolamide was used as a positive control. Relative carbonic anhydrase activity (%) = (catalytic rate of esterase reaction with inhibitor) / (catalytic rate of esterase reaction without inhibitor) × 100. IC₅₀ values of inhibitors were determined by plotting the percentage of enzyme activity against the inhibitor concentration (Bijari et al., 2015).

HPLC Conditions:

Chromatographic separation was performed on LC800 (GL Sciences). The auto-sampler and column heater temperature was maintained at 25°C and the injection volume was set at 15 µL for all analyses. The chromatographic separation was achieved on Agilent Zorbax Eclipse C18 column (50 × 2.1 mm, 1.7µm). The mobile phases consisted of acetonitrile (A) and water (B) both containing 0.1% formic acid. The gradient profile was set at 10% B from 0 to 1 min, 30% B at 8 min, 40% B at 12 min, 80% B at 16min, 95% B from 20-27 min, and finally 10% B at 28-35min. The flow rate for all separations was set at 0.7 mL/min.

MS/MS Conditions:

The sample was analyzed in hybrid TripleTOF® 5600 fitted with a DuoSpray™ ion source (AB Sciex, Concord, Canada). Every sample was injected twice in positive and negative polarity. The instrument was set to perform one TOF MS survey scan (150 ms) and 20 MS/MS scans (50 ms each) with a total duty cycle time of 1.2 s. The mass range of both scan types was 50-1000 m/z. Acquisition of MS/MS spectra was controlled by IDA function of the Analyst® TF software (AB Sciex, Concord, Canada).




Figure 1. Workflows for targeted and non targeted analysis

RESULTS

Effect of *L. siceraria* on carbonic anhydrase activity:

The LS metanolic extract (LSME) and three other fractions (aqueous, ethyl acetate, hexane) were assayed for carbonic anhydrase inhibition activity. The IC₅₀ values of aqueous fraction (LSAF) and acetazolamide (positive control) was estimated to be 403.34 ± 5.04 µg/mL, 220.32 ± 3.56 µg/mL respectively. (Fig. 1). It has been found that, LSAF inhibited carbonic anhydrase activity up to 89% at the concentration of 900 µg/ml in dose dependent manner.(Fig. 2). The result also suggested that LSAF inhibited carbonic anhydrase, reversibly. The enzyme kinetic analysis revealed mixed inhibition mode (i.e., both competitive and non-competitive) of LSAF, further supported by Lineweaver-Burk plot analysis of the kinetic data.




Fig. 2. IC₅₀ values of LS extract and fractions





Fig 3. Dose dependent inhibition of LSAF and ACZ

Identification of chemical compounds

#		Name	Formula	Isotope	Mass (Da)	Adduct	Extinction Mass (Da)	Width (Da)	Fragment Mass (Da)	Found A Mass (Da)	Error (ppm)	Expected RT (min)	RT Width (min)	Found A RT (min)	RT Delta (min)	Intensity	Threshold (cps)
✓	●●●●	1	C19H18O8	0	302.10017	-	301.0208	0.01	301.0208	301.0999	2.6	4.74	1	4.74	0	33895	
✓	●●●●	2	C19H18O8	0	300.08452	-	299.07724	0.01	299.07724	299.0777	1.9	8.22	1	8.22	0	15496	
✓	●●●●	3	C22H24O10	0	441.13015	-	441.12987	0.01	441.12987	441.1296	0.5	10.98	1	10.98	0	43481	
✓	●●●●	4	C22H24O10	0	442.14033	-	441.13676	0.01	441.13580	441.1303	0.8	12.13	1	12.13	0	15589	
✓	●●●●	5	C22H24O10	0	440.12030	-	440.11911	0.01	440.1182	440.1182	-0.7	10.96	1	10.96	0.01	6067	
✓	●●●●	6	C22H24O10	0	442.11921	-	441.10584	0.01	441.1058	441.1053	-0.8	13.19	1	13.19	0.08	1816	
✓	●●●●	7	C23H26O10	0	462.1758	-	461.16532	0.01	461.1653	461.164	-1.2	12.37	1	12.37	0	5977	
✓	●●●●	8	C23H26O10	0	461.17077	-	460.16430	0.01	460.16349	460.16349	-1.6	8.63	1	8.63	0	7663	
✓	●●●●	9	C37H40O17	0	755.2202	-	755.2197	0.01	755.2191	755.2191	-0.6	12.03	1	12.02	0.01	537	
✓	●●●●	10	C37H40O17	0	756.2272	-	755.22584	0.01	755.2256	755.2256	-4.6	12.01	1	11.99	0.02	154	
✓	●●●●	11	C10H10O4	0	193.04679	-	193.04679	0.01	193.04679	193.04679	0	15.14	1	15.14	0	2964	
✓	●●●●	12	C9H8O3	0	164.0474	-	163.04007	0.01	163.04004	163.04004	-4.7	12.15	1	12.15	0	27347	
✓	●●●●	13	C23H26O10	0	756.23212	-	755.22584	0.01	755.22546	755.22546	-4.6	12.01	1	11.99	0.02	154	
✓	●●●●	14	C23H26O10	0	460.0626	-	459.0498	0.01	459.03956	459.03956	-3.2	12.23	1	12.23	0.08	1930	
✓	●●●●	15	C23H26O10	0	462.06922	-	461.0687	0.01	461.06789	461.06789	-1.1	10.86	1	10.86	0	3385	
✓	●●●●	16	C23H26O10	0	461.07817	0.01	460.16349	0.01	460.16349	460.16349	-1.6	8.63	1	8.63	0	7663	
✓	●●●●	17	C19H17O4	0	314.0791	-	313.05003	0.01	313.05128	313.05128	3.4	15.14	1	15.14	0	2964	
✓	●●●●	18	C18H12O3	0	300.0786	-	299.07137	0.01	299.07121	299.07121	-4.1	14.8	1	14.8	0	46758	
✓	●●●●	19	C18H12O3	0	300.0786	-	301.04004	0.01	301.04004	301.04004	4.7	12.15	1	12.15	0	27347	
✓	●●●●	20	C18H12O3	0	301.07996	-	299.07069	0.01	299.07069	299.07069	-2.9	14.5	1	14.5	0	2670	
✓	●●●●	21	C9H8O3	0	162.01918	-	161.04002	0.01	161.04587	161.04587	133.3	22.98	1	23.27	0.29	74	