

# Quantitation of Pegylated Drug Conjugate Interferon Alfa-2b in Human Serum using QTRAP® 6500 System



Faraz Rashid<sup>1</sup>, Dipankar Malakar<sup>1</sup> Anoop Kumar<sup>1</sup> Manoj Pillai<sup>1</sup>, Praveen Kumar V<sup>2</sup>, Manoj Bob<sup>2</sup>, Ravisekhar K<sup>2</sup> <sup>1</sup>SCIEX, 121 DHR, Udyog Vihar, Phase IV, Gurgaon, Haryana, India; <sup>2</sup>Lupin Bioresearch Center, Pune Maharashtra, India

## INTRODUCTION

One of the most effective adaptable and acceptable procedure for stabilization of therapeutic protein is PEGylation which improves the shelf life. PEGylated therapeutics proteins are known for its stability and some of them are commercially available. There is a demand for developing sensitivities and selective LC-MS/MS based method for the quantitation which will help to understand the pharmacokinetics of these molecules. In this study we have developed MRM based quantitation method in human serum for Interferon Alfa-2b which is PEGylated therapeutic protein drug conjugate. In this method surrogate peptides were generated using trypsin digestion from the therapeutic protein. Various sample preparation and cleanup parameters were evaluated and optimized for achieving the highest sensitivity for PEGylated Interferon Alfa-2b (PEG-IFN)

## MATERIALS AND METHODS

### PEGylated Interferon Alfa-2 b (PEG-IFN)

PEGylated Interferon Alfa-2 b was provided by Lupin Bioresearch Center, Pune Maharashtra, India, at the concentration of 5ng/ul, subsequently diluted as per the study reported for quantitation.

### LC-MSMS Conditions:

#### LC gradient run for Information dependent Acquisition (IDA) to select surrogate peptide

A Shimadzu Nexra LC system with Xbridge dC18, 100x 2.1 mm, 3.5µm column at 40°C with an elution gradient of 5-35% acetonitrile (0.1% formic acid) with a 35 min gradient was used at a flow rate of 300µL/min. The injection volume was set to 10µL. Mobile Phase A : Water :01.% Formic Acid Mobile Phase B : ACN 01.% Formic Acid

### Gradient for MRM

A Shimadzu Nexra LC system with Xbridge dC18, 100x 2.1 mm, 3.5µm column at 45°C with a mobile phase gradient A (Water with 01.% FA ) and B (Methanol 0.1% FA) was used at a flow rate of 400µL/min. The injection volume was set to 30µL. Gradient is listed in Table No: 1

### Extraction of PEGylated proteins from Serum

Most of the proteins in serum samples were precipitated with water-miscible organic solvents and protein of interested (PEG-IFN) were collected in supernatant by centrifugation.

### LC-MS/MS Conditions:

An SCIEX QTRAP® 6500 LC/MS/MS system with Turbo V™ source and Electrospray Ionization (ESI) probe was used for quantitation. QTRAP based IDA workflow was used to identify a surrogate and sensitive peptide of PEGylated Interferon Alfa-2b and sequence confirmation of the same done by ProteinPilot™ 5.0 . Skyline software was used to identify the correct precursor and its fragments and the acquisition method was automatically created and exported to Analyst® 1.6 software for the respective peptide. Source and compound related parameters were optimized to achieve highest sensitivity. Three MRM transitions were selected for quantitation and identification of PEG-IFN in serum.

Time(min)	Flow (µl/min)	% A (0.1% FA in Water)	% B (0.1% FA in MeOH)
0.01	400	80	20
0.50	400	80	20
15.0	400	40	60
16.0	400	20	80
18.0	400	20	80
18.5	400	80	20
20.0	400	80	20


Table 1: LCMS gradient table for separation of Peptides in matrix used for MRM




SCIEX QTRAP® 6500



IonDrive™ Turbo V Source, IonDrive™ QJet Ion Guide, IonDrive™ High Energy Detector




### Information Dependent Acquisition (IDA) with XIC of selected peptide




IDA Method Criteria

Source Parameters		Compound	
CUR	30	DP	85
IS	5500	EP	10
TEM	500	CE	32
GS1	55	CXP	12
GS2	50	Dwell (ms)	400
CAD	06		



Source and compound Parameters for MRM




Schematic workflow for the Sample Preparation of Interferon-α2b in serum



Linearity of extracted PEG-IFN peptide in the range of 246 pg/ml-50,000 pg/ml in serum



Chromatograms of blank matrix and LLOQ



Calibration curve of peptide “SLSTNLQESLR” m/z 741.382+>1047.5

## CONCLUSION

- Highly specific and robust method is developed for PEGylated drug conjugate Interferon Alfa-2b.
- Doubly charge peptide were selected for better and robust ionization for MRM m/z 741.412+> 1047.54
- Liquid liquid Extraction (LLE) based method was optimized for the quantitation of PEG-IFN in serum matrix
- Linearity within a concentration range of 246-50,000 pg/ml of PEG-IFN was achieved in serum matrix for quantitation
- All signal responses were linear in the tested concentration range, represented by regression coefficients r<sup>2</sup> 0.992
- The limit of quantification (LOQ) was determined for target peptide in serum was estimated as the lowest concentration measured with % CV< 20%.

## REFERENCES

- Rinat R. et al., International Journal of Mass Spectrometry 312 (2012) 135– 143.
- Wu S, et al, Rapid Commun Mass Spectrom. 2011;25:281-90.
- Yang Z, et al., J Chromatogr B Analyt Technol Biomed Life Sci. 2009;877:1737-42.
- Long Yuan and Mingshe Zhu, Mater Methods 2015;5:133

## TRADEMARKS/LICENSING

For Research Use Only. Not for use in diagnostic procedures. AB Sciex is doing business as SCIEXThe trademarks mentioned herein are the property of Sciex Pte. Ltd. or their respective owners. AB SCIEX™ is being used under license. © 2017 Sciex. RUO-MKT-10-4183-A