

High-sensitivity detection of forensic drug panel in human whole blood

Using the SCIEX Triple Quad™ 7500 LC-MS/MS – QTRAP® Ready System, powered by SCIEX OS Software

Pierre Negri¹ and Ian Moore²

¹SCIEX, USA; ²SCIEX, CA

Achieving low levels of detection while maintaining reliable quantification is a key performance indicator to any robust bioanalytical workflow. In the forensic laboratory, the ability to accurately quantify large panels of analytes in complex biological matrices over a wide range of concentration is challenging and often requires dilution of the samples to bring them within the calibration range of the instrument. In order to further reduce sample preparation and re-analysis time, sufficient data sampling across the chromatographic peaks and over a wide linear dynamic range is required to ensure comprehensive analyte coverage and rugged quantitative performance.

In this technical note, the SCIEX Triple Quad 7500 System is presented as a robust quantitative platform with exceptional performance in speed, linear dynamic range and sensitivity for the accurate quantification of a panel of 49 drugs in human whole blood. The optimized method maximizes the linear dynamic range capabilities of the assay while retaining the required levels of accuracy and performance. The addition of new hardware features of the OptiFlow® Pro Ion Source, the D Jet™ Ion Guide and the E Lens™ Technology enhance ion sampling and desolvation,¹ resulting in unparalleled sensitivity and quantification limit improvements for the suite of drugs targeted in this study.



Key advantages of the SCIEX Triple Quad 7500 System for forensic drug panel analysis

- The SCIEX Triple Quad 7500 System is a robust instrument with unparalleled sensitivity and quantitative performance
- SCIEX OS Software provides an easy to use and intuitive platform for both for data acquisition and processing
- Fast (6.5 minutes) chromatographic run time achievable with the Scheduled MRM™ Algorithm Pro in SCIEX OS Software and fast polarity switching during data acquisition
- Robust and easy to use OptiFlow Pro Ion Source provides efficient sample ionization with intelligent source design eliminating the need for physical optimization
- Improved desolvation and ion sampling of system enhanced sensitivity, with LLOQ in the sub ng/mL to pg/mL range, demonstrated in complex biological matrices for the panel of drugs targeted in this study
- Accurate quantification of all the drugs in the panel, over a wide range of concentrations, without any sacrifice to data quality
- Overall performance of the system resulted in excellent precision and reproducibility in the reported data, even at low concentration levels




Figure 1. High linearity and linear dynamic range (LDR) demonstrated for chemically diverse panel of drugs. Calibration curves resulting from the calibration series of alprazolam (0.01-100 ng/mL) and codeine (0.01-50 ng/mL) showing excellent linear response even at pg/mL concentrations.

Experimental details

Target analytes and solutions: A total of 49 drugs and 18 deuterated internal standards were purchased from Cerilliant Corporation (Round Rock, TX). Two solutions were prepared in methanol: a standard mixture containing the 49 target analytes and an internal standard mixture containing the 18 deuterated internal standards. Table 1 lists the name, the calibration range, linear correlation value (R^2), and LLOQ, as well as the accuracy and precision reported at the LLOQ for each of the 49 target analytes used in this panel.

Calibrator preparation: Eight levels of calibrators ranging from 10 $\mu\text{g/mL}$ to 1 pg/mL were prepared in methanol. A 10 ng/mL IS standard stock solution was prepared in MeOH/water (20:80, v/v) for sample reconstitution prior to injection.

Sample preparation: 10 μL of each calibrator solution was spiked into 90 μL of human whole blood. Each spiked human whole blood sample was extracted by using a protein precipitation procedure. In short, 900 μL of methanol/acetonitrile (50:50, v/v) were added into each of the spiked human whole blood samples and vortexed for 1 min then followed by 3 min sonication and another 1 min of vortex mixing. The samples were then centrifuged for 5 min at 8,000 rpm. The supernatant was transferred out to a glass tube and completely dried down under nitrogen gas. The residues were reconstituted with 500 μL of a 10 ng/mL IS standard stock solution in methanol/water (20:80, v/v). The protein precipitation procedure is shown in Figure 2.

Liquid chromatography: HPLC separation was performed on an ExionLC™ System using a Phenomenex Kinetex Phenyl-Hexyl column (50 \times 2.1 mm, 2.6 μm , 00B-4495-E0). The separation conditions were identical to those previously described in a technical note.² Mobile phases were ammonium formate in water (MPA) and formic acid in methanol (MPB). The injection volume was 5 μL and the LC runtime was 6.5 min.

Mix	•10 μL of std mixture with 90 μL of human whole blood
Load to tube	•900 μL of MeOH: MeCN (50:50, v/v)
Vortex	•Vortex vigorously for 1 min
Sonicate	•Sonicate for 3 min
Vortex	•Vortex vigorously for 1 min
Centrifuge	•5 min at 8,000 rpm
Transfer	•Transfer supernatant to glass vial
Evaporate	•Evaporate to dryness under nitrogen
Reconstitute	•Add 500 μL of a 10 ng/mL IS std mixture in MeOH: water (20:80, v/v)

Figure 2. Protein precipitation procedure for human whole blood samples. A 9-step protein precipitation protocol was used for selectively extracting drugs from human whole blood samples for analysis with the SCIEX Triple Quad 7500 System.

Mass spectrometry: A SCIEX Triple Quad 7500 System equipped with an OptiFlow Pro Ion Source with an OptiFlow analytical probe and E Lens™ Technology was used. The ionization source was operated in electrospray ionization (ESI) mode in both positive and negative polarities. A single acquisition method consisting of 134 MRM transitions (98 for the drugs and 36 for the internal standards) was created using the Scheduled MRM Algorithm Pro in SCIEX OS Software 2.0. Two MRM transitions were monitored for each of the targeted drugs and each sample was injected in triplicate to build a data analysis processing method.

Data analysis: Data processing was performed using SCIEX OS Software. Detection and integration of the peaks from the background was achieved within the viewing window using the AutoPeak algorithm. Quantitative analysis was performed in the Analytics module of the software where calibration curves, concentration calculations, assay precision and accuracy statistics were generated.

Method development and optimization

A diluted, neat standard mixture containing the 49 target analytes was used for initial method development. The Scheduled MRM Algorithm Pro in SCIEX OS Software was used to optimize data sampling across each peak while maintaining optimal dwell times for each MRM transition to ensure reliable integration, quantification and confirmation of the peak for each target analyte. In addition, fast polarity switching was used to provide maximum analyte coverage. Most MRM transitions had 15 or more data points across each of the LC peaks, with 10 being the minimum number of data points across each peak for each of the 49 target analytes used in this study. Figure 3 shows the elution profile for the 49 targeted drugs resulting from the optimized data acquisition method.




Figure 3. Chromatographic profile of the 49 drugs targeted in this study. Extracted ion chromatograms (XICs) resulting from the optimized data acquisition using a neat standard mixture. Method optimization using the Scheduled MRM™ Algorithm Pro in SCIEX OS Software enabled collection of optimal data quality even during regions of the chromatogram when MRM concurrency was very high.

Optimized detection method enables accurate and reliable drug quantification

Control human whole blood samples spiked with the 49 target analytes were prepared at concentrations ranging from 1 µg/mL down to 0.1 pg/mL. Detection and integration of the peaks was performed automatically using the AutoPeak Algorithm in the Analytics module of SCIEX OS Software. Analyte concentration and ion ratio were calculated automatically in the software.

The potential for drastic variation in detected drug levels in toxicology case samples necessitates the use of an instrument with high range of linear response. Figure 1 shows calibration curves for two of the drugs targeted in this study. Concentration range of 4 (from 0.01-100 ng/mL) and 3.7 (from 0.01-50 ng/mL) orders linear dynamic range was demonstrated for alprazolam and codeine, respectively. Excellent linearity was observed across the concentration ranges analyzed with R^2 average values of 0.99761 and 0.99901 for alprazolam and codeine, respectively. Similar trends were observed for the other analytes used in this study.

OptiFlow Pro Ion Source and E Lens Technology leads to enhanced sensitivity

Developing robust workflows that can deliver high levels of sensitivity is critical to any toxicology laboratory needing to quantify a wide concentration of drugs extracted from biological specimens. To this end, the sensitivity of the SCIEX Triple Quad 7500 System was assessed by determining the lower limit of quantification (LLOQ) values for the 49 targeted drugs in the panel. LLOQ values were determined as the lowest concentration calibration level fitting the following standard performance requirements: signal-to-noise ratio (S/N) of at least 10, calculated concentration accuracy within 20% of 100%, precision (%bias) below 25%, and falling on a linear calibration curve with an R^2 value of at least 0.98.

Figure 4 shows the extracted ion chromatogram (XIC) traces and resulting calibration curves of the two MRM transitions monitored for 7-hydroxymitragynine (Figure 4A) and acetyl fentanyl (Figure 4B). The two series of XIC traces for both the quantifier and qualifier ions of each of the two detected drugs showed a high level of sensitivity and precision across the calibration series for concentrations ranging from 0.1 to 100 ng/mL for 7-hydroxymitragynine and 0.02-50 ng/mL for acetyl fentanyl, respectively. Six levels of calibrators were used to determine the ion ratio criteria for the quantifier and qualifier ions of these two targeted drugs. The results demonstrated excellent correlation of the generated regression curves covering concentration ranges far exceeding typical bioanalytical requirements.




Figure 4: High sensitivity and linearity for selected drugs in the forensic panel. A) Calibration curves and XIC traces resulting from the calibration series for the two transitions of 7-hydroxymitragynine from 0.1 to 100 ng/mL. B) Calibration curves and XIC traces resulting from the calibration series for the two transitions of acetyl fentanyl from 0.02 to 50 ng/mL. The calibration curves and XIC traces demonstrate excellent linearity and sensitivity.

Table 1 summarizes the statistical results obtained for the 49 analytes spiked in human whole blood. The table includes calibration range, linear correlation coefficient (R^2 Value), and LLOQ, as well as the accuracy and precision for each of the two MRM transitions monitored for each drug. Overall, the assay showed excellent reproducibility, precision, accuracy, and linearity, proving the robustness and performance of the developed method.

Figure 5 shows the statistical results and the calibration curves resulting from the peak area integration of dihydrocodeine (Figure 5A) and noroxycodone (Figure 5B) from 0.5 to 100 ng/mL. Excellent linearity, reproducibility, accuracy and precision was observed across the six calibration levels covering the concentration range. The assay showed excellent precision and accuracy, and the averaged R^2 values for the quantifier and qualifier ions were 0.99571 and 0.99323, respectively. Full quantification was achieved with SCIEX OS Software, designed for quick, intuitive and streamlined data processing with accurate and reliable results.




Figure 5. Statistical results and calibration curves for selected drugs in the forensic panel. Statistics pane and calibration curves for dihydrocodeine (A) and noroxycodone (B) from 0.5 to 100 ng/mL. Both analytes showed excellent linearity, reproducibility, accuracy and precision across the six calibration levels, proving the overall robustness of the method.

Conclusions

An optimized and sensitive method for the detection of a panel of 49 drugs in human whole blood is described using the SCIEX Triple Quad 7500 System. The use of the Scheduled MRM Algorithm Pro in SCIEX OS Software enabled optimization of data sampling. The addition of new hardware features, the OptiFlow Pro Ion Source, the D Jet Ion Guide and the E Lens, provided sensitive quantification of all the targeted drugs in the sub ng/mL range, with some down to the pg/mL levels, while maintaining linearity, precision and accuracy of measurement. This exceptional sensitivity was achieved without any sacrifice or compromise to data quality, as demonstrated by the excellent precision and accuracy observed at the LLOQ. Overall, the combination of the features on the SCIEX Triple Quad 7500 System results in unparalleled sensitivity improvement for the suite of drugs targeted in this study.

References

1. Enabling new levels of quantification - using the SCIEX Triple Quad™ 7500 LC-MS/MS System – QTRAP® Ready, powered by SCIEX OS Software. [SCIEX Technical Note RUO-MKT-02-11886-A](#).
2. High Sensitivity and Dynamic Range for 93-Compound Forensic Panel Analysis in Urine. [SCIEX Technical Note RUO-MKT-02-9914-A](#).



Table 1. Statistical results for the 49 drugs targeted in this workflow. The table includes calibration range, linear correlation coefficient (R^2 Value), and LLOQ, as well as the accuracy and precision at the LLOQ for each of the two MRM transitions monitored for each drug.

Compound	Calibration Range (ng/mL)	Linear Correlation (R^2)	LLOQ (ng/mL)	Accuracy at LLOQ (%)	Precision at LLOQ (%)
6-MAM 1	0.2-200	0.99191	0.2	82.62	2.74
6-MAM 2	0.2-200	0.99020	0.2	83.02	4.91
7-Aminoclonazepam 1	0.5-100	0.99910	0.5	117.72	2.66
7-Aminoclonazepam 2	0.5-100	0.99936	0.5	112.71	7.54
7-Hydroxymisragynine 1	0.1-100	0.98891	0.1	84.76	1.72
7-Hydroxymisragynine 2	0.1-100	0.99006	0.1	90.58	4.33
Acetyl fentanyl 1	0.02-50	0.98272	0.02	113.43	6.33
Acetyl fentanyl 2	0.02-50	0.98393	0.02	112.01	9.49
Alpha-Hydroxylprazolam 1	1-100	0.99918	1	93.38	2.57
Alpha-Hydroxylprazolam 2	1-100	0.99817	1	93.67	7.67
Alprazolam 1	0.01-100	0.99780	0.01	84.15	5.96
Alprazolam 2	0.01-100	0.99741	0.01	80.98	1.33
Amphetamine 1	0.2-200	0.98925	0.2	117.13	22.79
Amphetamine 2	0.2-200	0.98784	0.2	86.02	8.95
Benzoyllecgonine 1	0.5-100	0.9777	0.5	100.00	5.82
Benzoyllecgonine 2	0.5-100	0.99776	0.5	103.70	3.07
Buprenorphine 1	0.4-400	0.99018	0.4	98.61	24.21
Buprenorphine 2	0.4-400	0.98312	0.4	96.09	5.49
Carisoprodol 1	2-200	0.95267	2	95.22	13.59
Carisoprodol 2	2-200	0.98049	2	117.40	12.85
Codeine 1	0.01-50	0.99897	0.01	80.29	13.48
Codeine 2	0.01-50	0.99904	0.01	93.47	17.80
Dextromethorphan 1	0.5-100	0.99598	0.5	119.17	10.95
Dextromethorphan 2	0.5-100	0.99772	0.5	101.42	5.86
Diazepam 1	0.5-100	0.99916	0.5	118.51	1.07
Diazepam 2	0.5-100	0.99522	0.5	118.30	0.96
Dihydrocodeine 1	0.5-100	0.98889	0.5	98.40	5.04
Dihydrocodeine 2	0.5-100	0.98006	0.5	117.31	4.45
EDDP 1	0.1-200	0.98973	0.1	90.93	8.74
EDDP 2	0.1-200	0.99128	0.1	109.63	13.35
Fentanyl 1	0.4-400	0.98210	0.4	112.74	3.66
Fentanyl 2	0.4-400	0.98615	0.4	110.68	2.30



Compound	Calibration Range (ng/mL)	Linear Correlation (R ²)	LLOQ (ng/mL)	Accuracy at LLOQ (%)	Precision at LLOQ (%)
Gabapentin 1	2-200	0.98159	2	115.14	5.53
Gabapentin 2	2-200	0.98326	2	113.60	3.94
Hydrocodone 1	1-100	0.97652	1	99.30	5.59
Hydrocodone 2	1-100	0.99070	1	103.96	7.82
Hydromorphone 1	1-100	0.98541	1	110.69	0.35
Hydromorphone 2	1-100	0.98618	1	97.29	1.25
Lorazepam 1	0.5-100	0.99699	0.5	0.6.76	2.89
Lorazepam 2	0.5-100	0.99680	0.5	119.49	5.71
MDA 1	10-100	0.92279	10	109.46	19.09
MDA 2	10-100	0.99617	10	108.54	10.94
MDEA 1	2-200	0.99567	2	100.00	7.62
MDEA 2	2-200	0.99682	2	100.00	5.95
MDMA 1	5-200	0.9855	5	108.39	1.94
MDMA 2	5-200	0.99562	5	102.96	5.53
Methadone 1	1-100	0.99040	1	97.53	4.64
Methadone 2	1-100	0.99330	1	115.54	5.37
Methamphetamine 1	2-200	0.99361	2	112.10	324
Methamphetamine 2	2-200	0.99550	2	105.70	1.71
Methylphenidate 1	1-100	0.99276	1	96.82	0.65
Methylphenidate 2	1-100	0.98423	1	87.52	4.34
Midazolam 1	1-100	0.99861	1	92.81	3.91
Midazolam 2	1-100	0.99484	1	114.07	3.21
Mitragynine 1	0.2-200	0.99572	0.2	97.85	16.06
Mitragynine 2	0.2-200	0.99422	0.2	96.69	6.71
Morphine 1	1-100	0.98897	0.2	116.82	0.73
Morphine 2	1-100	0.98830	0.2	116.50	1.46
Naloxone 1	1-100	0.98323	1	100.71	0.90
Naloxone 2	1-100	0.98563	1	97.88	4.11
Naltrexone 1	0.5-100	0.99273	0.5	119.45	4.17
Naltrexone 2	0.5-100	0.98868	0.5	119.72	18.77
Norbuprenorphine 1	0.4-200	0.97745	0.4	98.35	4.73
Norbuprenorphine 2	0.4-200	0.97920	0.4	98.05	0.50
Norcodeine 1	0.5-100	0.98892	0.5	119.46	2.73
Norcodeine 2	0.5-100	0.99030	0.5	117.30	1.03
Nordiazepam 1	1-100	0.99073	1	107.85	0.48
Nordiazepam 2	1-100	0.99087	1	107.46	1.11

Compound	Calibration Range (ng/mL)	Linear Correlation (R ²)	LLOQ (ng/mL)	Accuracy at LLOQ (%)	Precision at LLOQ (%)
Norfentanyl 1	0.4-100	0.99366	0.4	100.28	4.03
Norfentanyl 2	0.4-100	0.99601	0.4	99.91	4.13
Norhydrocodone 1	1-100	0.98777	1	93.83	6.16
Norhydrocodone 2	1-100	0.99486	1	94.06	7.91
Noroxycodone 1	0.5-100	0.99411	0.5	110.52	5.61
Noroxycodone 2	0.5-100	0.99235	0.5	112.35	5.62
Norpropoxyphene 1	1-200	0.97989	1	105.35	9.60
Norpropoxyphene 2	1-200	0.98323	1	104.65	10.70
O-Desmethyltramadol 1	1-100	0.96310	1	93.34	13.91
O-Desmethyltramadol 2	1-100	0.96164	1	90.77	4.82
Oxazepam 1	0.5-100	0.99704	0.5	85.31	7.74
Oxazepam 2	0.5-100	0.99195	0.5	86.37	14.15
Oxycodone 1	0.5-100	0.98486	0.5	113.73	11.84
Oxycodone 2	0.5-100	0.98492	0.5	115.65	10.38
Oxymorphone 1	1-100	0.98986	0.5	114.44	2.78
Oxymorphone 2	1-100	0.99246	0.5	116.20	1.98
PCP 1	0.5-100	0.99167	0.5	103.36	10.99
PCP 2	0-100	0.98672	0.5	92.84	14.21
Pregabalin 1	1-200	0.98380	1	94.53	16.48
Pregabalin 2	1-200	0.98277	1	94.79	16.12
Tapentadol 1	0.5-100	0.98018	0.5	94.01	13.36
Tapentadol 2	0.5-100	0.99810	0.5	103.41	7.72
Temazepam 1	0.05-100	0.99217	0.05	113.97	4.89
Temazepam 2	0.05-100	0.99249	0.05	119.15	4.38
Tramadol 1	0.5-100	0.99674	0.5	102.19	10.33
Tramadol 2	0.5-100	0.99779	0.5	110.82	3.63
Zolpidem 1	1-100	0.98577	1	93.87	4.83
Zolpidem 2	1-100	0.99810	1	103.41	7.81
THC-COOH 1	1-100	0.9537	1	116.92	5.86
THC-COOH 2	1-100	0.9846	1	118.37	7.68

The SCIEX clinical diagnostic portfolio is For In Vitro Diagnostic Use. Rx Only. Product(s) not available in all countries. For information on availability, please contact your local sales representative or refer to <https://sciex.com/diagnostics>. All other products are For Research Use Only. Not for use in Diagnostic Procedures.

Trademarks and/or registered trademarks mentioned herein, including associated logos, are the property of AB Sciex Pte. Ltd. or their respective owners in the United States and/or certain other countries.

© 2020 DH Tech. Dev. Pte. Ltd. RUO-MKT-02-11427-A. AB SCIEX™ is being used under license.



Headquarters
500 Old Connecticut Path | Framingham, MA 01701 USA
Phone 508-383-7700
sciex.com

International Sales
For our office locations please call the division
headquarters or refer to our website at
sciex.com/offices